

→ PROTEOR QUATTRO 2™ PATIENT TRIAL / DELIVERY CHECKLIST



GENERAL INFORMATION

Date: _____

Facility: _____

CP: _____

Patient ID: _____

Pyramid or Threaded Top: ☐ P ☐ T

Foot: _____ Side: ☐ R ☐ L Size: _____ Cat: _____ Toe: ☐ Sandal ☐ Reg Footshell Color: ☐ L ☐ M ☐ D ☐ B

FACILITY CHECKLIST

- Access to Ramps or Slopes
- Access to Stairs
- Access to Uneven Terrain

FOR DELIVERY START AT STEP 3

STEP 1: SET EXPECTATIONS (PATIENT AND CP)

- The QUATTRO 2 is a revolutionary design and will not walk like your current knee
- Discuss how the valve control technology is constantly making changes to optimize gait on every step or situation
- Consider taking a video of each knee to provide a comparison (Consent form)

STEP 2: ACTIVITIES ON CURRENT KNEE

Prior to the patient trialing the QUATTRO 2, have the patient perform these tasks on their current knee:

- ☐ Walking Backwards
- ☐ Side steps
- ☐ Slow, Medium, Fast walking
- ☐ Transition from slow to fast, fast to slow
- ☐ Stairs and Ramps step over step (if they currently do)
- ☐ Sitting and sit to stand
- ☐ Uneven Terrain (if available)

STEP 3: INSTALL AND ALIGN THE QUATTRO 2

- ☐ Ensure proper alignment prior to programming the knee
- ☐ Bench Alignment (no Stance Flexion: 0-5mm anterior; with Stance Flexion 0-5mm posterior to knee center)
- ☐ Static Alignment- patient standing quietly (no Stance Flexion: 0-10mm anterior; with Stance Flexion 0-10mm posterior to knee center)
- ☐ Dynamic Alignment- patient walking

STEP 4: COMPLETE PROGRAMMING IN GAITLAB

- ☐ Calibration (consider walking outside the parallel bars with walking aid for Step 3)
- ☐ Resistance Settings
- ☐ Demonstrations

→ PROTEOR QUATTRO 2™ PATIENT TRIAL / DELIVERY CHECKLIST

STEP 5: DEMONSTRATIONS HIGHLIGHTING FEATURES

- ☐ Walking Backwards
- ☐ Side steps
- ☐ Slow, Medium, Fast walking
- ☐ Transition from slow to fast, fast to slow
- ☐ Stairs and Ramps step over step (if they currently do)
- ☐ Sitting and sit to stand
- ☐ Uneven Terrain (if available)

	KEYPAD	GAITLAB	FI APP
Stance Lock (on/off if applicable)	—	☐	☐
Obstacle Assist	—	☐	—
Flexion Lock	☐	☐	☐
Sitting Mode	—	☐	☐
Select Activities	—	☐	☐
Setting Alternate Activities	—	☐	☐
Switch Activities (Walking and Alternate)	☐	—	—
Add New Activities	—	☐	☐
Adjust Activities	—	☐	☐
Cycling	—	☐	☐
Golfing	—	☐	☐
Sprinting	—	☐	☐
Elliptical	—	☐	☐
Rehabilitation	—	☐	☐
Remote Firmware Update	—	☐	☐
Cadence and Usage Reports	—	☐	☐
Outcome Measures	—	☐	☐

FINAL SETTINGS

Sitting: _____ Heel Rise: _____ Swing Ext: _____ Stance Ext: _____ Stair/Ramp: _____

Obstacle Assist: ☐ On ☐ Off Sitting Mode: ☐ On ☐ Off Stance Lock: ☐ On ☐ Off

Activities: ☐ Cycling ☐ Golfing ☐ Sprinting ☐ Elliptical ☐ Rehabilitation ☐ Custom: _____

Notes: _____